

Looking for help?

Find answers to your questions

Search..



🔒 Xtrodes	5
🔒 Technical & Research Information	37
🔒 Troubleshooting and Cleaning	16
🔒 Headsets	11
🔒 Software	14

📄 Connecting to Service issue (<https://www.emotiv.com/tools/knowledge-base/software/connecting-to-service-issue>)

📄 Where do I find EMOTIV Software to download? (<https://www.emotiv.com/tools/knowledge-base/software/where-do-i-find-emotiv-software-to-download>)

▼ EmotivPRO

📄 Emotiv data for qEEG interpretation (<https://www.emotiv.com/tools/knowledge-base/emotivpro/emotiv-data-for-qeeg-interpretation>)

📄 Why is my EEG quality flickering in EmotivPRO?
(<https://www.emotiv.com/tools/knowledge-base/emotivpro/why-is-my-eeeg-quality-flickering-in-emotivpro>)

📄 How to submit Feedback via EmotivPRO for issue investigation
(<https://www.emotiv.com/tools/knowledge-base/emotivpro/how-to-submit-feedback-via-emotivpro-for-issue-investigation>)

📄 What is EEG Quality in the EmotivPRO? (<https://www.emotiv.com/tools/knowledge-base/emotivpro/what-is-eeeg-quality-in-the-emotivpro>)



- Does EmotivPRO automatically remove artifacts from the EEG data collected?
(<https://www.emotiv.com/tools/knowledge-base/emotivpro/does-emotivpro-automatically-remove-artifacts-from-the-eeeg-data-collected>)
- How can timing markers be entered into the data stream?
(<https://www.emotiv.com/tools/knowledge-base/emotivpro/how-can-timing-markers-be-entered-into-the-data-stream>)

EMOTIV® ([/](#)) [Q \(/search\)](#) [🛒 \(/shop.emotiv.com/\)](#)

➤ BrainVIZ

🔒 FAQ	25
🔒 Licenses	3
🔒 Xtrodes	24
🔒 Account Settings	2
🔒 Billing and Payment	2
🔒 Shipping, Taxes, and Delivery	13
🔒 Subscription	7
🔒 Refunds, Returns and Replacements	7

Home (<https://www.emotiv.com>) ➤ Docs (<https://www.emotiv.com/tools/knowledge-base>) ➤
EmotivPRO (<https://www.emotiv.com/tools/knowledge-base/category/emotivpro>) ➤
Understanding Sensor Values in Exported CSV Data
(<https://www.emotiv.com/tools/knowledge-base/emotivpro/understanding-sensor-values-in-exported-csv-data>)

Understanding Sensor Values in Exported CSV Data

Please refer to the EmotivPRO User Manual (<https://emotiv.gitbook.io/emotivpro-v2-0/>) for a detailed explanation. The values are expressed in microvolts but converted directly from a 14-bit unsigned integer. The DC level is offset by approximately 4100 uV, but the absolute value of the offset can be affected by slow DC drift so we do not attempt to correct it. Please apply a high-pass filter to the data – we recommend a characteristic cutoff around 0.16 Hz or greater. You can refer to this article (<https://emotiv.gitbook.io/emotivpro-v3/notes-on-the-data/dc-offset>) for how to remove the DC offset from the data.

As a side note, you can also use our cloud computing analyzer (<https://www.emotiv.com/emotivpro/analyze/>) to effortlessly post-process data to analyze your data. On the Analyzer page, you can select basic EEG processing or Basic EEG processing and Transformation.

"Basic EEG Processing" applies a slew rate limit, a high-pass filter, and a common-mode noise removal step (subtract signals which are common to a set of channels = background body potential signals that managed to pass through the analog filtering systems). Users can choose parameters at each step. This returns cleaner EEG signals with the DC offset removed. If you choose "Basic EEG Processing and Transformation", Analyzer performs the above filtering, then applies FFT and returns average power in each specified frequency band.

Updated on 10 Jul 2025

What are your Feelings?



Product

Solutions

Support


Company








© 2025 EMOTIV, All rights reserved. (https://www.shopify.com?utm_campaign=poweredby&utm_medium=shopify&utm_source=onlinestore)

[Do Not Sell My Personal Information](#)

(<https://www.emotiv.com/pages/ccpa-opt-out>)

 Your Privacy Choices (Cookie Settings)

 (<https://www.facebook.com/emotiv>)  (<https://www.instagram.com/emotiv/>)  (<https://www.linkedin.com/company/emotiv>)  (<https://www.youtube.com/c/emotivstation>)  (<https://twitter.com/emotiv>)

*Disclaimer – EMOTIV products are intended to be used for research applications and personal use only. Our products are not sold as Medical Devices as defined in EU directive 93/42/EEC. Our products are not designed or intended to be used for diagnosis or treatment of disease.